RACTIVEONE

RECOVERY BOOTS · 4 CHAMBERS

REVIVE ACTIVE

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EN
User Manual

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INHALT / CONTENTS / SOMMAIRE

DE Deutsch

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
01 IN THE BOX

WELCOME TO THE TEAM


Revive Active brings compression recovery home. Four air chambers per leg build pressure from the bottom up – for fresh legs after hard sessions. The control unit sits right on the boot – everything within reach.

1

2

3

4


IN THE BOX

1	Boots	× 1 pair
2	Air hose	× 1
3	USB charging cable	× 1
4	Control unit	× 1

02 CONTROL PANEL

ALL ON THE BOOT

The control unit display shows program, pressure, time, battery and the chamber currently active.


1. **Pressure display**
 30–240 mmHg
2. **Program**
 Active program 1–5
3. **Power / Pause**
 On, off, pause
4. **MODE**
 Switch program
5. **Pressure –**
 Lower pressure
6. **Timer**
 15 / 30 / 45 min
7. **Pressure +**
 Raise pressure
8. **Battery**
 Charge level
9. **Time display**
 Remaining minutes
10. **Position**
 Active chamber 1–4
11. **Charging port**
 USB, 5 V ⎓ 1 A

03 MOUNTING

MOUNTING

Before you start: fit the control unit and connect both boots with the hose.

CONNECT CONTROL UNIT & HOSE

1. Slide the control unit into the holder on the boot from above until it clicks.
2. Push the hose plug firmly into the socket underneath the control unit – arrow facing up.
3. Connect the other end of the hose to the second boot.

04 QUICK START

GO IN 4 STEPS

GET READY

- Charge the control unit with the USB cable before first use.
- Connect the hose to the control unit and boots (chapter 03).
- Sit comfortably with legs extended and zip up the boots.

1 **Switch on**
Press the power button on the control unit.

2 **Set pressure**
Use – / + to choose between 30 and 240 mmHg.

3 **Choose time**
Press the timer button: 15 → 30 → 45 min.

4 **Choose program**
Press MODE: program 1 → 2 → 3 → 4 → 5.

TIP Start your first session at low pressure, e.g. 80 mmHg, and increase until it feels right.

05 PROGRAMS

YOUR RHYTHM

Five programs, each with its own sequence. Press MODE to cycle through them. Areas marked in blue are active.

1 **Wave**

One chamber after another fills – a wave from foot to thigh.

2 **Sequential**

Chambers fill one by one and hold pressure until the whole leg is wrapped.

3 **Double Wave**

Two chambers work together as a double wave moving up.

4 **All**

All chambers fill at once and release together.

5 **Sequential + Double Wave**

Building compression first, followed by a double wave.

PRESSURE LEVELS

30 · 50 · 80 · 120 · 160 · 200 · 240 mmHg


TIMER

15 · 30 · 45 min

06 AFTER TRAINING

PACK UP RIGHT

To store the boots, pull the hose plug off the control unit. Always press the clips on both sides of the plug – so it does not break.

Press the side clips, then pull the plug down.

CLEANING & CARE

- Switch the unit off before cleaning.
- Do not wash the boots – wipe with a clean, slightly damp cloth.
- Clean the display and control unit with a soft, dry cloth.
- Store clean and dry, away from children, heat and direct sunlight.

TROUBLESHOOTING

Unit does not start
Battery empty → **Charge with the USB cable, then switch on.**

Unit switches off
Set time has run out → **Simply restart.**

Boots barely inflate
Hose not plugged in properly or air leak → **Check the plug; if it persists, contact support.**

07 SAFETY & SPECS

BEFORE YOU START

- ! Do not unzip the boots while the unit is running – always switch off first.
- ! Only charge with the supplied USB cable.
- ! Do not open or repair the unit yourself.
- ! Do not use in wet areas or if liquid has entered the unit.
- ! Do not walk around wearing the boots during use.
- ! Keep away from sharp objects.
- ! If you feel unwell, stop immediately.
- ! Consult your physician before first use. This product is not intended to diagnose, treat or cure any disease.

DO NOT USE WITH

Acute venous thrombosis	Deep thrombophlebitis
Acute skin inflammation	Open wounds, burns
Heart disease	Unstable hypertension
Pacemaker	Pregnancy

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© 2026 Ractiveone