RACTIVEONE

RECOVERY BOOTS · 8 CHAMBERS

REVIVE ACTIVE PRO

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EN
User Manual

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INHALT / CONTENTS / SOMMAIRE

DE Deutsch

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01 IN THE BOX

WELCOME TO THE TEAM

Revive Active Pro brings professional compression recovery home. Eight air chambers per leg build pressure from the bottom up – so you show up fresher for your next session. Give this manual five minutes and you are ready to go.

Control unit with touch display


IN THE BOX

1	Control unit	× 1
2	Y-connector (1 to 2)	× 1
3	Air hose	× 2
4	Charger	× 1
5	Leg cuff	× 2
6	Carry case	× 1

02 CONTROL PANEL

FULL CONTROL

The touch display shows program, time and pressure at a glance. Adjust every setting with the buttons on the left (–) and right (+).

1. **Mode – / +**
 Choose program
2. **Time – / +**
 Set duration
3. **Pressure – / +**
 Set pressure
4. **Position**
 Select chamber
5. **Confirm**
 Confirm selection
6. **Battery**
 Shows charge level
7. **On / Pause**
 Start and pause


03 SETUP


PUT ON & CONNECT

1 Connect the air hose to both cuffs.

2 Push the hose plug firmly into the pump – arrow facing up.

3 Sit comfortably with legs extended and zip up the cuffs.

Insert the hose plug into the pump socket with the arrow facing up.

TIP Switch the pump on before you connect the cuffs.

04 QUICK START

9 STEPS TO YOUR SESSION

BEFORE FIRST USE

Press the power button. If the pump does not start, the battery is empty – plug in the charger. The LED on the adapter lights red and turns green while charging. When fully charged, the battery indicator shows all bars.

1 Hose to cuffs
Connect firmly until secure.

2 Hose to pump
Push the main plug in firmly.

3 Put on the cuffs
Fit comfortably and close the zip.

4 Switch on
Front switch to "I" – the pump goes to standby. "O" turns it off.

5 Set time: 10–90 min
Use Time – / +. Default 30 min. Hold the button for 5-minute steps. The pump switches off automatically when time is up.

6 Set pressure: 80–260 mmHg
Use Pressure – / +. Default 80 mmHg. Hold for faster changes. The pump stops below 80 mmHg.


7 Choose program
Use Mode – / +. Default is program 1. You can only switch while paused.

8 Skip chambers (optional)
Deactivate single chambers – see chapter 05.

9 Press On / Pause
Once everything is set, your session starts.

05 PROGRAMS

YOUR RHYTHM

Six programs, six sequences: the chart shows the order in which the eight chambers fill. Highlighted segments are active.

Switching programs? Press On / Pause first, then choose with Mode – / +.

SKIP SPECIFIC CHAMBERS

1. Press Position: the selection moves one chamber, the selected chamber flashes.
2. Press Confirm: the LED goes off – this chamber will not inflate.
3. Repeat to skip more chambers.


06 AFTER TRAINING

PACK UP RIGHT


To pack up, simply pull the hose plug out of the pump. Press the clips on both sides of the plug so it does not break.

Keep hose and cuffs connected – next time just plug the other end back into the pump. Frequent disconnecting wears out pins and hose.

When travelling you can detach the hose from the cuffs to save space. Always grip the pin – never pull the hose.

Press the side clips, then pull the plug down.

Grip the pin and pull up – not the hose.

CLEANING & CARE

- Switch the pump off before cleaning.
- Do not wash the pump or cuffs.
- Wipe with a dry or slightly damp cloth.
- Store clean and dry.
- Keep away from direct sunlight and heat sources (e.g. gas stove).

07 SAFETY & SPECS

BEFORE YOU START

- ! Do not unzip the cuffs while the unit is running – always switch off first.
- ! Do not get the pump wet.
- ! Place the pump on a flat surface.
- ! Do not walk around wearing the cuffs during use.
- ! Do not repair or service the unit while in use.
- ! Consult your physician before first use. This product is not intended to diagnose, treat or cure any disease.

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© 2026 Ractiveone