RACTIVEONE

HEAT & COLD MASSAGE GUN

REVIVE PULSE PRO

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EN
User Manual

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INHALT / CONTENTS / SOMMAIRE

DE Deutsch

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
01 OVERVIEW

WELCOME TO THE TEAM

Revive Pulse Pro releases tension deep in your muscles and fascia – before or after training. Five intensity levels plus heat and cold via the thermo head: loosen up faster and show up fresher for your next session.

A Control panel	**B** Massage head
C Battery indicator	**D** USB-C charging port

02 ATTACHMENTS

FIVE ATTACHMENTS

Hot & cold head
Gentle on large areas like lower legs, thighs and arms.

U-shaped head
For neck, back, arms and Achilles tendon.

Cylinder head
For palms, soles and deep muscle areas.


Soft head
For shoulders, back and thighs – extra gentle.

Ball head
For shoulders, wrists, thighs and lower legs.


CHANGING HEADS

1. Heads push straight in and pull straight out – no screwing.
2. On the hot & cold head, line up the red bar with the gap in the piston.
3. Always push the head in all the way.

03 OPERATION

FULL CONTROL

- **A** Vibration level indicator
- **B** On/off & pause
- **C** Heat/cold indicator
- **D** Heat button
- **E** Vibration button
- **F** Cold button

1 Switch on

Short press the power button: the panel lights up for 1 s and the light at the bottom of the handle turns on – the gun is in standby.

2 Start vibration

Short press the vibration button. 5 levels – each press moves up one level, the indicator shows the current one.

3 Pause

While running, short press the power button to pause vibration; press again to resume.

4 Switch off

Hold the power button for 2 seconds – works in any mode.

04 HEAT & COLD

HOT OR COLD

Heat and cold only work with the hot & cold head. You can also use vibration, heat and cold separately.

HEAT

In standby, press the heat button – indicator lights red. Press again = next level.

Level 1	**42 °C**
Level 2	**46 °C**
Level 3	**50 °C**

COLD

In standby, press the cold button – indicator lights blue. Press again = next level.

Level 1	**20 °C**
Level 2	**15 °C**
Level 3	**10 °C**

Depending on your surroundings, values may vary by up to ±2 °C (heat) or ±5 °C (cold).

AUTO STOP AFTER 15 MINUTES


Vibration, heat and cold each run for a maximum of 15 minutes – then Revive Pulse Pro switches off automatically. Want another round? Just restart.

SPEED PER LEVEL

LEVEL	HOT & COLD HEAD	OTHER HEADS
Level 1	**1600 RPM**	**1800 RPM**
Level 2	**1800 RPM**	**2100 RPM**
Level 3	**2000 RPM**	**2400 RPM**
Level 4	**2200 RPM**	**2700 RPM**
Level 5	**2400 RPM**	**3000 RPM**

05 CHARGING

FULLY CHARGED

1. Charge for approx. 5 hours before first use. When full, the charging indicator lights white and goes off after 5 minutes.
2. Plug the included USB-C cable into the port at the bottom of the handle. Any 5 V USB adapter works (recommended: 5 V / 2 A).
3. You can top up anytime – whatever the charge level.
4. Do not use while charging. Unplug the cable after charging.

BATTERY INDICATOR ON THE HANDLE

White	Battery full
Green	Battery medium
Red	Battery low – please charge
Red flashing	Almost empty – switches off after 5 flashes

06 HELP & CARE

TROUBLESHOOTING

Won't start	Charge the battery. If unused for a long time, charge for more than 8 hours.
Won't start in heat	Overheat protection above 40 °C ambient: let it rest below 40 °C for 2 hours, then restart.
Won't charge	Plug the cable in firmly. If device or charger is too hot, let them rest below 40 °C for 2 hours.
Battery deeply drained	Keep switched off, charge fully, unplug – then use as normal.
Head slips out	Push the head in all the way (see chapter 02).
Stops suddenly	Auto stop after 15 min, or protection because the head was blocked. Short press the power button and restart.

CLEANING

- Wipe with a slightly damp cloth.
- Dry with a soft, dry cloth.
- Never immerse in water.

STORAGE

- For long breaks, store clean in the original box.
- Keep dry, well ventilated and away from corrosive gases.
- When stored, charge once a month.
- Do not dispose of with household waste – take it to a recycling point.

07 SAFETY & SPECS

BEFORE YOU START

- ! Not for children under 8. Minors only after consulting a doctor.
- ! Use only on dry, clean skin or over clothing. Press gently and move for approx. 60 s per area.
- ! Do not immerse; keep water and objects out of the vents.
- ! Keep away from fire and heat sources. Do not drop or modify.
- ! Do not operate or charge unattended. Do not use while charging.
- ! Pinch hazard: keep fingers and hair away from the shaft and back of the head when switching on.
- ! Take care if your sensitivity to heat or cold is reduced.
- ! With muscle injuries, consult your doctor before use.

ASK YOUR DOCTOR FIRST IF YOU HAVE

Pregnancy, cancer, heart disease, a pacemaker or defibrillator, recent surgery, aneurysm, bleeding or blood thinners, neuropathy – and on sensitive areas such as the head and joints.

SPECIFICATIONS

Model	**B330**	Weight	**approx. 776 g**
Size	**202 × 52 × 169 mm**	Power	**35 W (max. 70 W)**
Input	**5 V ⎓ 2 A, USB-C**	Battery	**11.1 V / 2400 mAh**
Levels	**5**	Charge time	**approx. 4–6 h**
Run time vibration	**3–4 h**	Run time heat / cold	**8–9 h / 3–4 h**

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© 2026 Ractiveone