RACTIVEONE

HEAT & COOL RECOVERY

THERMOFLEX

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EN
User Manual

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INHALT / CONTENTS / SOMMAIRE

DE Deutsch

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
01 IN THE BOX

WELCOME TO THE TEAM

Thermoflex puts heat and cold for your recovery in your bag. Cool down after training, warm up to loosen – flexibly on neck, shoulders, knees or calves. Cordless, battery-powered and on in seconds.


IN THE BOX

1	Thermoflex device	× 1
2	Strap frame	× 1
3	Velcro straps	× 3
4	USB-C charging cable	× 1
5	User manual	× 1

02 DEVICE OVERVIEW

AT A GLANCE

1. **Cool button**
 Cooling on / change level
2. **Heat button**
 Heating on / change level
3. **USB-C charging port**
 Charges the battery
4. **Heat/cool surface**
 Sits directly on your skin
5. **Air inlet**
 Active in cooling mode
6. **Air outlet**
 Active in cooling mode

IMPORTANT IN COOLING MODE
Always keep the air inlet (5) and air outlet (6) clear – otherwise cooling performance drops.

03 USE & WEARING

WHERE TO USE IT

Head

Neck

Back

Thighs

Shoulders

Arms

Knees

Calves

HOW TO PUT IT ON

1 Attach the strap to one side of the device with the velcro.

2 Bend the device and place it on the target area.

3 Wrap the strap around and fix it on the other side.

TIP On neck, shoulders and knees, just bend the device to fit and put it in place.

04 STRAP FRAME & CHARGING

STRAP FRAME

On arms and legs, the strap frame gives Thermoflex an even more secure fit by changing where the strap is fixed.

1 Fit the strap frame on the device and attach the strap to one side.


2 Place the device and wrap the strap around.

3 Fix the velcro on the other side.

CHARGING

1 Fully charge before first use or after a long break: plug the USB-C cable into the charging port.

2 Connect the cable to a power adapter and plug it in.

3 The indicator flashes while charging.

4 Fully charged after approx. 3.5 h. Hold any button for 3 seconds – when the indicator goes off, unplug. Do not use while charging.

05 OPERATION

HOT OR COLD

HEAT

With the device off, hold the heat button for 3 s. Starts at level 1, red LED. Short press = change level.

Level 1	**38 ± 3 °C**
Level 2	**40 ± 3 °C**
Level 3	**42 ± 3 °C**

COOL

With the device off, hold the cool button for 3 s. Starts at level 1, blue LED. Short press = change level.

Level 1	**15 ± 3 °C**
Level 2	**10 ± 3 °C**
Level 3	**5 ± 3 °C**

The number of lit LEDs shows your current level: 1, 2 or 3.

HOT & COLD CYCLE MODE

With the device off, hold both buttons for 3 s – the LED flashes blue. The device switches between cool and heat every 3 minutes, starting with cool level 2. A short press on the matching button changes the level of the current phase.

SWITCHING OFF

Thermoflex switches off automatically after 15 minutes in every mode. Manually: hold any button for 1 second. Want to keep going? Just restart.

Cooling values are measured at 20 °C room temperature. They may vary slightly with your surroundings – that is normal.

06 HELP & CARE

TROUBLESHOOTING

Does not start	Check the battery and charge if needed. Hold the button for at least 1 second.
Stops during use	The device switches off after 15 min – just restart. Battery empty? Charge fully, then continue.
Won't charge / short run time	Only charge at 10–40 °C. Check that USB-C cable and adapter are firmly connected. Fully charge first time or after a long break. The device does not charge while running – switch off first.
Red skin	Normal after heat/cold – fades within 30 min. Check whether you react to the material.
Warm after switching off	Normal after cooling – returns to normal in a few minutes.

CARE

- Unplug the charging cable before cleaning.
- Do not use detergents, thinners or petrol.
- Wipe with a soft, dry cloth – do not rinse under water.
- Clean the contact surface after each use, with a little medical disinfectant if needed.

STORAGE

- Unplug the charging cable.
- Store dry and out of direct sunlight.
- For long breaks, charge fully and keep in a dry, ventilated place.
- Keep away from children and pets.
- Do not place on fabric, paper or delicate wood.

07 SAFETY & SPECS

BEFORE YOU START

- ! Do not use near water or fire.
- ! Do not use in strong electromagnetic fields, high heat or humidity.
- ! Do not use near flammable or explosive substances (e.g. petrol stations).
- ! Do not drop or open – contact our support if it is faulty.
- ! Do not plug or unplug the cable with wet hands.
- ! Do not use with other devices or more than 2 Thermoflex at once.
- ! Do not use while driving.
- ! Stop immediately if you feel unwell.
- ! Do not give to children or people who cannot express discomfort.
- ! Take care if you are sensitive to temperature changes, heat or cold.
- ! Contains a lithium battery: do not disassemble, modify or throw into fire. Dispose of according to local regulations.

DO NOT USE WITH

Pregnancy, high fever, malignant tumours, untreated heart conditions, acute illness, impaired sensation. Thermoflex is not a medical device and does not replace medical treatment.

SPECIFICATIONS

Model	**KB-CL801Pro**	Power	**20 W**
Size	**149 × 76 × 50 mm**	Charging voltage	**5 V**
Battery	**3200 mAh**	Port	**USB-C**
Weight	**387 g**	Charge time	**approx. 3.5 h**

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© 2026 Ractiveone