RACTIVEONE

HOT · COLD · RED LIGHT · VIBRATION

VIBRAGRIP

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EN
User Manual

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INHALT / CONTENTS / SOMMAIRE

DE Deutsch

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

01 IN THE BOX

WELCOME TO THE TEAM

VibraGrip combines deep vibration, heat, cold and red light in one device. Ideal for activating tight muscles before training or releasing them afterwards.

1


2

IN THE BOX

1	VibraGrip	× 1
2	USB-C charging cable	× 1

02 DEVICE OVERVIEW

YOUR VIBRAGRIP


1 Hot & cold flat head	2 LED display
3 Control panel	4 Ambient light
5 Air outlet	6 Air inlet

03 HEADS & PORTS

THREE HEADS, ONE DEVICE

Depending on what you need, use a different side of the head – flat, targeted or with red light.


7	Silicone massage head	8	Red light LEDs
9	Finger pressure head	10	Red light head
11	On/off button	12	USB-C charging port

04 CONTROL PANEL

EVERYTHING IN HAND

Display and buttons sit right on the handle. The display shows vibration level, heat, cold and battery.

1. Vibration indicator
2. Cooling indicator
3. Heating indicator
4. Battery indicator
5. Vibration +
6. Vibration –
7. Heat
8. Cold
9. Red light

05 HOW TO USE

HOW IT WORKS

1 Switch on

Briefly press the on/off button. The ambient light comes on – the device is ready. Without input it switches off after 1 minute.

2 Choose a function

Use the control panel buttons to set vibration, heat, cold or red light.

3 Switch off

Hold the on/off button for 2 seconds. The device switches off automatically after 15 minutes – simply restart.

VIBRATION

0–10

Adjust level by level with + / –.

HEAT

40 · 45 · 50 °C

Press briefly to switch level.

COLD

15 · 10 · 5 °C

Press briefly to switch level.

RED LIGHT

660 + 850 nm

Red and infrared light – press briefly for on/off.


NOTE

- Temperature depends on the environment and may vary by ±3 °C.
- In cold mode, keep the air inlet and outlet uncovered – otherwise cooling is reduced.

06 CHARGING

START FULLY CHARGED

Fully charge VibraGrip before first use or after a long break.

1 Plug the USB-C cable into the charging port on the bottom.

2 Connect the cable to a USB power adapter and plug it in.

3 While charging, the battery indicator on the display flashes.

4 After approx. 3 hours the indicator shows 5 bars – unplug the device.

NOTE Do not use the device while charging.

07 CARE & TROUBLESHOOTING

KEEP IT IN SHAPE

- Unplug the charging cable before cleaning.
- Wipe with a soft, dry cloth only – do not clean under running water.
- Do not use detergents, thinners, petrol or other chemicals.
- Do not open, repair or modify the device.
- Store dry, away from moisture, direct sunlight, children and pets.
- If not used for a long time, store fully charged in a dry, ventilated place.

TROUBLESHOOTING

Device does not start
Check the battery and charge. Press the on/off button only briefly.

No heat / vibration
Place the heat head flat on the skin. Check the battery.

Device stops during use
It switches off automatically after 15 minutes – restart. Charge if the battery is low.

Won't charge / short battery life
Charge at 10–40 °C. Check that cable and adapter are firmly connected. Fully charge before first use.

If the problem persists, disconnect the power and contact our support.

08 SAFETY & SPECS

BEFORE YOU START

- ! Read this manual completely before first use.
- ! Always unplug the charging cable before use and cleaning.
- ! Do not immerse in water or other liquids. If the device falls into water while charging: unplug immediately, do not touch the device.
- ! Contains a lithium battery – do not open, modify or throw into fire.
- ! Do not plug or unplug the charging cable with wet hands.
- ! Use the heat function with care, especially on sensitive skin. Do not use on the same spot for long – risk of burns.
- ! Children and people with limited abilities only under supervision. Do not leave within reach of unattended children.
- ! If you feel unwell, stop immediately.
- ! Not a substitute for medication or medical devices. Consult your physician if you have health concerns.

SPECIFICATIONS

Model	**HXR-B510**	Size	**216 × 78 × 91 mm**
Weight	**600 g**	Power	**20 W**
Battery	**7.4 V · 2200 mAh**	Charging	**5 V ⎓ 2 A**
Port	**USB-C**	Charge time	**approx. 3 h**

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© 2026 Ractiveone